



# Mini-spectrometer

[ TM series ]

C10082MD

C10083MD

**For UV to near IR, integrating optical system,  
image sensor and circuit**

TM series mini-spectrometers are polychromators integrated with optical elements, an image sensor and a driver circuit. Two models are available: C10082MD (TM-UV/VIS-MOS) and C10083MD (TM-VIS/NIR-MOS). Light to be measured is guided into the entrance port of TM series through an optical fiber and the spectrum measured with the built-in image sensor is output from the USB port to a PC for data acquisition. No external power supply is required since USB bus power is used for circuit operation.

Mini-spectrometer TM series comes supplied with free evaluation software that allows setting measurement conditions, acquiring and saving data, and displaying graphs. DLL is also supplied as accessory items to allow the users to configure their own measurement software.

## Features

- High throughput due to transmission grating made of quartz
- Highly accurate optical characteristics
- No external power supply required: uses USB bus power
- Wide spectral response range
- Easy to install into equipment
- Wavelength conversion factor\*1 is recorded in internal memory.
- Supports external trigger input\*2

\*1: A conversion factor for converting the image sensor pixel number into a wavelength is recorded in the module. A calculation factor for converting the A/D converted count into the input light intensity is not provided.

\*2: Coaxial cable for external trigger input sold separately. Refer to the "Mini-spectrometers Selection Guide" for details on external triggers.

## Applications

- Light source spectrum measurement
- Sunlight or illumination analysis
- Absorbance measurement

## Optical characteristics

| Parameter  | TM-UV/VIS-MOS | TM-VIS/NIR-MOS | Unit  |
|--|---------------|----------------|-------|
|  | C10082MD      | C10083MD       |       |
| Spectral response range                              | 200 to 800    | 320 to 1000    | nm    |
| Spectral resolution (Spectral response half width)*3 | 6 max.        | 8 max.         | nm    |
| Wavelength reproducibility*4                         | -0.2 to +0.2  |                | nm    |
| Wavelength temperature dependence                    | -0.4 to +0.4  |                | nm/°C |
| Spectral stray light*3 *5                            | -35 max.      | -33 max.       | dB    |

\*3: Depends on the slit opening. Values were measured with the slit listed in the table "Structure / Absolute maximum ratings".

\*4: Measured under constant light input conditions

\*5: When monochromatic light of the following wavelengths is input, spectral stray light is defined as the ratio of the count measured at the input wavelength, to the count measured in a region of the input wavelength  $\pm 40$  nm.

C10082MD: 500 nm, C10083MD: 650 nm

## Electrical characteristics

| Parameter                         | Specification | Unit |
|-----------------------------------|---------------|------|
| A/D conversion                    | 16            | bit  |
| Integration time                  | 5 to 10000    | ms   |
| Interface                         | USB 1.1       | -    |
| USB bus power current consumption | 100 max.      | mA   |

## Structure / Absolute maximum ratings

| Parameter                           | Value                                  | Unit   |
|-------------------------------------|--|--------|
| Dimensions (W × D × H)              | 94 × 90 × 55                           | mm     |
| Weight                              | 470                                    | g      |
| Image sensor                        | CMOS linear image sensor (S8378-1024Q) | -      |
| Number of pixels                    | 1024                                   | pixels |
| Slit* <sup>6</sup> (H × V)          | 70 × 800                               | μm     |
| NA* <sup>7</sup>                    | 0.22                                   | -      |
| Connector for optical fiber         | SMA905D                                | -      |
| Operating temperature* <sup>8</sup> | +5 to +40                              | °C     |
| Storage temperature* <sup>8</sup>   | -20 to +70                             | °C     |


\*6: Entrance slit aperture size

\*7: Numerical aperture (solid angle)

\*8: No condensation


Note: Exceeding the absolute maximum ratings even momentarily may cause a drop in product quality. Always be sure to use the product within the absolute maximum ratings.

## Spectral resolution vs. wavelength (typical example)



KACC00265EA


## Linearity (typical example)



Note: A/D output is the output with dark output is subtracted when light is input. The difference between the ideal value and typical example contains a measurement error. The smaller the A/D output, the larger the measurement error.

KACC00244EA

### Dark output vs. ambient temperature (typical example)




Note: A/D output is the sum of the sensor and circuit offset outputs and the sensor dark output.

KACCB0245EA

### Optical component layout


TM series mini-spectrometers use a transmission holographic grating made of quartz and precision optical components arranged on a rugged optical base, making it possible to deliver high throughput and highly accurate optical characteristics.



KACCC0287EA

### Connection example (transmission light measurement)

Light to be measured is guided into the entrance port of TM series through an optical fiber and the spectrum measured with the built-in image sensor is output through the USB port to a PC for data acquisition. There are no moving parts inside the unit so stable measurements are obtained at all times. An optical fiber that guides light input from external sources allows a flexible measurement setup.



KACCC0288EF


### ❖ Evaluation software package (supplied with unit)

Installing the evaluation software package (Spec Evaluation. exe)\*9 into your PC allows running the following basic tasks:

- Measurement data acquisition and save
- Measurement condition setup
- Module information acquisition  
(wavelength conversion factor, polychromator type, etc.)
- Graphic display
- Arithmetic operation
- Pixel number to wavelength conversion
- Comparison calculation with reference data (transmittance, reflectance)
- Dark subtraction
- Gaussian approximation (peak position and count, FWHM)

Note:

- Two or more mini-spectrometers can be connected and used with one PC simultaneously.
- The external trigger input function works with DLL, but does not function on the supplied evaluation software. If using an external trigger input, the user software must be configured to support that function.



\*9: Compatible OS: Microsoft® Windows® XP Professional SP3 (32-bit)\*10  
 Microsoft® Windows® Vista Business SP2 (32-bit)\*10  
 Microsoft® Windows® 7 Ultimate SP1 (32-bit)\*10  
 Microsoft® Windows® 7 Ultimate SP1 (64-bit)\*10

DLL for controlling hardware is also provided.

You can develop your own measurement programs by using a following software development environment.


Microsoft® Visual Studio® 2008 (SP1) Visual C++®\*10

Microsoft® Visual Studio® 2008 (SP1) Visual Basic®\*10

\*10: Microsoft, Windows, Visual Studio, Visual C++ and Visual Basic are either registered trademarks or trademarks of Microsoft Corporation in the United States and other countries.


### ❖ Measurement example (C10082MA)

Line spectra from visible LED



KACCB0126EB

### ▣ Dimensional outline (unit: mm, tolerance unless otherwise noted: $\pm 0.5$ )



Weight: 470 g

KACCA0171EE

### ▣ Accessories

- USB cable
- Dedicated software (evaluation software, sample software, DLL)

### ▣ Options (sold separately)

- Coaxial cable for external trigger input A10670
- Optical fibers for light input

| Type no. | Product name                                     | Core diameter ( $\mu\text{m}$ ) | Specification   |
|----------|--|---------------------------------|---|
| A9762-01 | Fiber for UV/visible range<br>(resistance to UV) | 600                             | NA=0.22, length 1.5 m,<br>connectorized SMA905D<br>at both ends |



The C10082MD and C10083MD conform to the European EMC directives  
(Applied standard: EN 61326-1 Class B).

## Mini-spectrometer lineup

| Type no.                | Type                           |  | Spectral response range (nm)  |            |     |     |      |      |      |      |      |      |      |      |      |  | Spectral resolution<br>max.<br>(nm) | Image sensor |                         |   |
|-------------------------|--------------------------------|--|-------------------------------|------------|-----|-----|------|------|------|------|------|------|------|------|------|--|-------------------------------------|--------------|-------------------------|---|
|                         |                                |  | 200                           | 400        | 600 | 800 | 1000 | 1200 | 1400 | 1600 | 1800 | 2000 | 2200 | 2400 | 2600 |  |                                     |              |                         |   |
| C10082CA                | TM series                      | TM-UV/VIS-CCD<br>High sensitivity            |                               |            |     |     |      |      |      |      |      |      |      |      |      |  |                                     |              | 6                       | Back-thinned type<br>CCD image sensor           |
| C10082CAH               |                                | TM-UV/VIS-CCD<br>High resolution             |                               | 200 to 800 |     |     |      |      |      |      |      |      |      |      |      |  |                                     | 1*           |                         |   |
| C10082MD                |                                | TM-UV/VIS-MOS<br>Wide dynamic range          |                               |            |     |     |      |      |      |      |      |      |      |      |      |  |                                     |              | 6                       | CMOS linear<br>image sensor                     |
| C10083CA                |                                | TM-VIS/NIR-CCD<br>High sensitivity           |                               |            |     |     |      |      |      |      |      |      |      |      |      |  |                                     |              | 8<br>(λ=320 to 900 nm)  | Back-thinned type<br>CCD image sensor           |
| C10083CAH               |                                | TM-VIS/NIR-CCD<br>High resolution            |                               |            |     |     |      |      |      |      |      |      |      |      |      |  |                                     |              | 1*<br>(λ=320 to 900 nm) |   |
| C10083MD                |                                | TM-VIS/NIR-MOS<br>Wide dynamic range         |                               |            |     |     |      |      |      |      |      |      |      |      |      |  |                                     |              | 8                       | CMOS linear<br>image sensor                     |
| <div>NEW</div> C11697MA |                                | TM-VIS/NIR-MOS-II<br>Trigger-compatible      |                               |            |     |     |      |      |      |      |      |      |      |      |      |  |                                     |              | 8                       | CMOS image sensor<br>with amp array             |
| C9404CA                 |                                | TG series                                    | TG-UV-CCD<br>High sensitivity |            |     |     |      |      |      |      |      |      |      |      |      |  |                                     |              |                         | 3   |
| C9404CAH                | TG-UV-CCD<br>High resolution   |  |                               |            |     |     |      |      |      |      |      |      |      |      |      |  |                                     |              | 1*                      | Back-thinned type<br>CCD image sensor           |
| <div>NEW</div> C9405CB  | TG-SWIR-CCD-II<br>IR-enhanced  |  |                               |            |     |     |      |      |      |      |      |      |      |      |      |  |                                     |              | 5<br>(λ=550 to 900 nm)  | IR-enhanced<br>back-thinned CCD<br>image sensor |
| <div>NEW</div> C11713CA | TG-RAMAN-I<br>High resolution  |  |                               |            |     |     |      |      |      |      |      |      |      |      |      |  |                                     |              | 0.3*                    | Back-thinned type<br>CCD image sensor           |
| <div>NEW</div> C11714CA | TG-RAMAN-II<br>High resolution |  |                               |            |     |     |      |      |      |      |      |      |      |      |      |  |                                     |              | 0.3*                    | Back-thinned type<br>CCD image sensor           |
| C9406GC                 | TG series                      | TG-NIR<br>Non-cooled type                    |                               |            |     |     |      |      |      |      |      |      |      |      |      |  |                                     |              | 7                       | InGaAs linear<br>image sensor                   |
| C9913GC                 |                                | TG-cooled NIR-I<br>Low noise (cooled type)   |                               |            |     |     |      |      |      |      |      |      |      |      |      |  |                                     |              | 7                       |   |
| C9914GB                 |                                | TG-cooled NIR-II<br>Low noise (cooled type)  |                               |            |     |     |      |      |      |      |      |      |      |      |      |  |                                     |              | 8                       |   |
| C11118GA                |                                | TG-cooled NIR-III<br>Low noise (cooled type) |                               |            |     |     |      |      |      |      |      |      |      |      |      |  |                                     |              | 20                      |   |
| C11007MA                | RC series                      | RC-VIS-MOS<br>Spectrometer module            |                               |            |     |     |      |      |      |      |      |      |      |      |      |  |                                     |              | 9                       | CMOS linear<br>image sensor                     |
| C11008MA                |                                | RC-SWIR-MOS<br>Spectrometer module           |                               |            |     |     |      |      |      |      |      |      |      |      |      |  |                                     |              | 8                       | IR-enhanced<br>CMOS linear<br>image sensor      |

\* Typ.

## For installation into mobile measuring equipment

| Type no. | Type      |                                  | Spectral response range (nm) |     |     |     |      |      |      |      |      |      |      |      |      |   | Spectral resolution max. (nm)        | Image sensor |
|----------|-----------|----------------------------------|------------------------------|-----|-----|-----|------|------|------|------|------|------|------|------|------|---|--------------------------------------|--------------|
|          |           |                                  | 200                          | 400 | 600 | 800 | 1000 | 1200 | 1400 | 1600 | 1800 | 2000 | 2200 | 2400 | 2600 |   |                                      |              |
| C11009MA | RC series | RC-VIS-MOS<br>Spectrometer head  |                              |     |     |     |      |      |      |      |      |      |      |      |      | 9 | CMOS linear image sensor             |              |
| C11010MA |           | RC-SWIR-MOS<br>Spectrometer head |                              |     |     |     |      |      |      |      |      |      |      |      |      | 8 | IR-enhanced CMOS linear image sensor |              |

## Ultra-compact type for installation into mobile measuring equipment

| Type no.                | Type      |                                    | Spectral response range (nm) |     |     |     |      |      |      |      |      |      |      |      |      |  | Spectral resolution<br>max.<br>(nm) | Image sensor                |
|-------------------------|-----------|------------------------------------|------------------------------|-----|-----|-----|------|------|------|------|------|------|------|------|------|--|-------------------------------------|-----------------------------|
|                         |           |                                    | 200                          | 400 | 600 | 800 | 1000 | 1200 | 1400 | 1600 | 1800 | 2000 | 2200 | 2400 | 2600 |  |                                     |                             |
| C10988MA                | MS series | MS-VIS-MOS<br>Spectrometer head    |                              |     |     |     |      |      |      |      |      |      |      |      |      |  | 14                                  | CMOS linear<br>image sensor |
| <div>NEW</div> C11708MA |           | MS-SWNI R-MOS<br>Spectrometer head |                              |     |     |     |      |      |      |      |      |      |      |      |      |  | 20                                  |                             |

Information described in this material is current as of September, 2012.

Product specifications are subject to change without prior notice due to improvements or other reasons. Before assembly into final products, please contact us for the delivery specification sheet to check the latest information.

Type numbers of products listed in the delivery specification sheets or supplied as samples may have a suffix "(X)" which means preliminary specifications or a suffix "(Z)" which means developmental specifications.

The product warranty is valid for one year after delivery and is limited to product repair or replacement for defects discovered and reported to us within that one year period. However, even if within the warranty period we accept absolutely no liability for any loss caused by natural disasters or improper product use.

Copying or reprinting the contents described in this material in whole or in part is prohibited without our prior permission.

**HAMAMATSU**

www.hamamatsu.com

HAMAMATSU PHOTONICS K.K., Solid State Division

1126-1 Ichino-cho, Higashi-ku, Hamamatsu City, 435-8558 Japan, Telephone: (81) 53-434-3311, Fax: (81) 53-434-5184

U.S.A.: Hamamatsu Corporation, 360 Foothill Road, P.O.Box 6910, Bridgewater, N.J. 08807-0910, U.S.A., Telephone: (1) 908-231-0960, Fax: (1) 908-231-1218

Germany: Hamamatsu Photonics Deutschland GmbH, Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49) 8152-375-0, Fax: (49) 8152-265-8

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trappu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: 33-(1) 69 53 71 00, Fax: 33-(1) 69 53 71 10

United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire AL7 1BW, United Kingdom, Telephone: (44) 1707-294888, Fax: (44) 1707-325777

North Europe: Hamamatsu Photonics Norden AB: Thorshamnsgatan 35 16440 Kista, Sweden, Telephone: (46) 8-509-031-00, Fax: (46) 8-509-031-01

Italy: Hamamatsu Photonics Italia S.R.L.: Strada della Moia, 1 int. 6, 20020 Arese, (Milano), Italy, Telephone: (39) 02-935-81-733, Fax: (39) 02-935-81-741

China: Hamamatsu Photonics (China) Co., Ltd.: 1201 Tower B, Jiaming Center, No.27 Dongsanhuan Beilu, Chaoyang District, Beijing 100020, China, Telephone: (86) 10-6586-6006, Fax: (86) 10-6586-2866