



### **Mini-spectrometer**

TG series

C9406GC, C9913GC, C9914GB

# For near IR, integrating optical system, image sensor and circuit

HAMAMATSU TG series mini-spectrometers are polychromators integrated with optical elements and an image sensor. Light to be measured is guided into the entrance port of TG series through an optical fiber and the spectrum measured with the built-in image sensor is output from the USB port to a PC for data acquisition. Non-cooled type and cooled type are provided. Non-cooled type is a palmtop-size unit and operates on USB bus power. Cooled type allows accurate measurement with low noise by cooling the image sensor. Two models are available: C9913GC (TG-cooled NIR-I) and C9914GB (TG-cooled NIR-II). The TG series comes supplied with evaluation software that allows setting measurement conditions, acquiring and saving data, and displaying graphs. DLL is also supplied as accessory items to allow the users to configure their own measurement software.

#### Features

- High throughput due to transmission grating made of quartz
- Highly accurate optical characteristics
- G9406GC: No external power supply required (uses USB bus power) \*1
- Low noise measurement (Cooled type: C9913GC, C9914GB)
- Compact design for easy assembly
- Wavelength conversion factor\*2 is recorded in internal memory

#### Applications

C9406GC (TG-NIR)

- Water content measurement
- **■** Optical communication component testing
- Film thickness measurement

C9913GC (TG-cooled NIR-I), C9914GB (TG-cooled NIR-II)

- Water content measurement
- Component analysis in food, agriculture fields, etc.
- Process control for chemical products
- \*1: C9913GC, C9914GB: Each requires 5 V and 12 V power supplies.
- \*2: A conversion factor for converting the image sensor pixel number into a wavelength is recorded in the module. A calculation factor for converting the A/D converted count into the input light level is not provided.

#### Optical characteristics

Parameter	TG-NIR	TG-cooled NIR-I	TG-cooled NIR-II	Unit					
Parameter	C9406GC	C9406GC							
Spectral response range	900 to 1700	900 to 1700	1100 to 2200	nm					
Spectral resolution (FWHM)*3	7 max.	8 max.	nm						
Wavelength reproducibility*4	-0.2 to +0.2	-0.2 to +0.2	-0.4 to +0.4	nm					
Wavelength temperature dependence	-0.02 to +0.02	-0.02 to +0.02	-0.04 to +0.04	nm/°C					
Spectral stray light*3 *5	ral stray light* <sup>3</sup> * <sup>5</sup> -35 max.								

<sup>\*3:</sup> Depends on the slit opening. Values were measured with the slit listed in the table "-Structure / Absolute maximum ratings".

<sup>\*4:</sup> Measured under constant light input conditions

<sup>\*5:</sup> When monochromatic light of the following wavelengths is input, spectral stray light is defined as the ratio of the count measured at the input wavelength, to the count measured in a region of the input wavelength ±40 nm. C9406GC, C9913GC: 1300 nm, C9914GB: 1650 nm

#### Electrical characteristics

Parameter	C9406GC	C9913GC	C9914GB	Unit
A/D conversion		bit		
Integration time	5 to 1	ms		
Interface		USB 1.1		-
Current consumption of USB bus power		mA		


#### **Structure / Absolute maximum ratings**


Parameter	C9406GC	C9914GB	Unit					
Dimensions (W $\times$ D $\times$ H)	38.5 × 106 × 86	mm						
Weight	270	1700	1700	g				
Image sensor	InGaAs linear image sensor (G9204-512D)	TE-cooled type InGaAs linear image sensor (G9204-512S)	TE-cooled type InGaAs linear image sensor	-				
Number of pixels*6	512	512	256	pixels				
Slit*7 (H × V)	70 × 500	70 × 500	70 × 250	μm				
NA*8		-						
Connector for optical fiber			-					
Operating temperature*9	+5 to +40	+5 to +35 (-	+5 to +30*10)	°C				
Storage temperature*9		-20 to +70		°C				
Power supply for cooling element max.*11	-	5/2.8	V/A					
Power supply for cooling fan	-	12/0.2						

Note: Exceeding the absolute maximum ratings even momentarily may cause a drop in product quality. Always be sure to use the product within the absolute maximum ratings.


- \*6: No defective pixel (at low gain). Defective pixels are those whose electrical and optical characteristics do not meet our specifications.
- \*7: Entrance slit aperture size
- \*8: Numeric aperture (solid angle)
- \*9: No condensation
- \*10: For controllable cooling temperature
- \*11: Maximum value in steady state. Note that inrush current flows at start-up. The connector for connection to cooling element and cooling fan power supply is attached (C9913GC, C9914GB).


#### Spectral resolution vs. wavelength (typical example)






#### Linearity (typical example)






Note: A/D output is the output with dark output is subtracted when light is input. The difference between the ideal value and typical example contains a measurement error. The smaller the A/D output, the larger the measurement error.




Note: A/D output is the output with dark output is subtracted when light is input. The difference between the ideal value and typical example contains a measurement error. The smaller the A/D output, the larger the measurement error.

#### Dark output vs. ambient temperature (typical example)




Note: A/D output is the sum of the sensor and circuit offset outputs and the sensor dark output.



Note: A/D output is the sum of the sensor and circuit offset outputs and the sensor dark output.

#### - Optical component layout


TG series mini-spectrometers use a transmission holographic grating made of quartz and precision optical components arranged on a rugged optical base, making it possible to deliver high throughput and highly accurate optical characteristics.



KACCC0256EA

#### Connection example (transmission light measurement)

Light to be measured is guided into the entrance port of TG series through an optical fiber and the spectrum measured with the built-in image sensor is output through the USB port to a PC for data acquisition. There are no moving parts inside the unit so stable measurements are obtained at all times. An optical fiber that guides light input from external sources allows a flexible measurement setup.



\* External power supply should be prepared by the user. C9406GC: No external power supply required (uses USB bus power)

#### Evaluation software package (supplied with unit)


Installing the evaluation software package (Spec Evaluation.exe)\*12 into your PC allows running the following basic tasks:

- · Measurement data acquisition and save
- · Measurement condition setup
- · Module information acquisition (wavelength conversion factor, polychromator type, etc.)
- · Graphic display
- Arithmetic operation

Pixel number to wavelength conversion Comparison calculation with reference data (transmittance, reflectance) Dark subtraction

Gaussian approximation

(peak position and count, FWHM)



#### Note:

- · Two or more mini-spectrometers can be connected and used with one PC simultaneously.
- The external trigger input function does not work with the evaluation software. If using an external trigger input or designing original application software, the user software must be configured to support that function.
- \*12: Compatible OS: Microsoft® Windows® XP Professional SP3 (32-bit)\*13

Microsoft® Windows® Vista Business SP1 (32-bit)\*13

Microsoft® Windows® 7 Ultimate SP1 (32-bit)\*13

Microsoft® Windows® 7 Ultimate SP1 (64-bit)\*13

DLL for controlling hardware is also provided.

You can develop your own measurement programs by using a following software development environment.

Microsoft® Visual Studio® 2008 (SP1) Visual C++®\*13

Microsoft® Visual Studio® 2008 (SP1) Visual Basic®\*13


\*13: Microsoft, Windows, Visual Studio, Visual C++ and Visual Basic are either registerd trademarks or trademarks of Microsoft Corporation in the United States and other countries.



#### Measurement example

Film thickness measurement (white light interferometry)

Thickness of 10  $\mu m$  thick food wrapping film (polyvinylidene chloride) was measured with C9406GC (TG-NIR).



#### Note:

Principle of film thickness measurement:

In film thickness measurement utilizing white light interferometry, an interference spectrum resulting from internal reflections between the front and back surfaces of a film is obtained.

The film thickness can then be determined by calculation from the spectral peak count, wavelength range, refractive index of film and incident light angle.

## Dimensional outline (unit: mm, tolerance unless otherwise noted: ±0.5)

C9406GC (2×) M3 tap depth 5.0 from backside 85.6 85.6 49.7 20.0 HAMAMATSU 16.0 41.8 68.6 38.5 105.5 Weight: 270 g

KACCB0095EA

KACCA0146ED

#### **Mini-spectrometer**

TG series C9406GC, C9913GC, C9914GB

#### C9913GC, C9914GB 218 142 ◈ COOLED NIR-II (λ: 2.2) Spectrometer COOLED NIR-II (λ: 2.2) Spectrometer 80 � ◈ ♦

Weight: 1.7 kg

KACCA0158EC

#### Accessories

- · USB cable
- · Dedicated software (evaluation software, sample software, DLL)
- · External power supply connector [made by LEMO S.A.: FGG0B304CLAD56 (only available in C9913GC and C9914GB)]

#### Options (sold separately)

Optical fibers for light input

Type no.	Product name	Specification
A9763-01		Core diameter 600 µm, NA=0.22, length 1.5 m,
A3703-01	Optical fiber for visible/fical infrared range	connectorized SMA905D at both ends

### CE

The C9406GC, C9913GC, and C9914GB conform to the European EMC directives. (applied standard: EN 61326-1 Class B)

#### Mini-spectrometer lineup

Type no		Tune							Spe	ctral i	es	ponse	rar	nge	(nr	n)							Spectral resolution	T
Type no.		Туре	200	) 4	00	600	8	300	100	120	00	1400	16	00	180	0Ó	2000	) 2	200	24	100	2600	max. (nm)	Image sensor
C10082CA		TM-UV/VIS-CCD High sensitivity																					6	Back-thinned type
C10082CAH		TM-UV/VIS-CCD High resolution		20	00 t	o 80	0																1*	CCD image sensor
C10082MD	Se	TM-UV/VIS-MOS Wide dynamic range																					6	CMOS linear image sensor
C10083CA	1 series	TM-VIS/NIR-CCD High sensitivity																					8 (λ=320 to 900 nm)	Back-thinned type
C10083CAH	Σ	TM-VIS/NIR-CCD High resolution			2	20 to	10	200															1* (λ=320 to 900 nm)	CCD image sensor
C10083MD		TM-VIS/NIR-MOS Wide dynamic range			٥.	20 10	) 10																8	CMOS linear image sensor
NEW C11697MA		TM-VIS/NIR-MOS-II Trigger-compatible																					8	CMOS image sensor with amp array
C9404CA		TG-UV-CCD High sensitivity	21	10 to 401	n																		3	Back-thinned type CCD image sensor
C9404CAH	SE	TG-UV-CCD High resolution	20	יטיד טו טיסי	U																		1*	Back-thinned type CCD image sensor
NEW C9405CB	series	TG-SWNIR-CCD-II IR-enhanced				50	00 t	o 1	100														5 (λ=550 to 900 nm)	IR-enhanced back-thinned CCD image sensor
NEW C11713CA	2	TG-RAMAN-I High resolution					500	to	600														0.3*	Back-thinned type CCD image sensor
NEW C11714CA		TG-RAMAN-II High resolution							79	0 to	92 	0											0.3*	Back-thinned type CCD image sensor
C9406GC		TG-NIR Non-cooled type								900	to	1700											7	
C9913GC	series	TG-cooled NIR-I Low noise (cooled type)										1700											7	InGaAs linear
C9914GB	TG s	TG-cooled NIR-II Low noise (cooled type)										11	.00	to	220	00							8	image sensor
C11118GA		TG-cooled NIR-III Low noise (cooled type)											90	0 to	o 2	55(	)						20	
C11007MA	series	RC-VIS-MOS Spectrometer module			340	to 7	80																9	CMOS linear image sensor
C11008MA	RC s	RC-SWNIR-MOS Spectrometer module					640	to	1050														8	IR-enhanced CMOS linear image sensor

<sup>\*</sup> Typ.

For installation into	For installation into mobile measuring equipment																
Type no.		Туре	200	400	600	800			ponse 1400			2000	2200	2400	2600	Spectral resolution max. (nm)	Image sensor
C11009MA	eries	RC-VIS-MOS Spectrometer head		340	to 78	30										9	CMOS linear image sensor
C11010MA	RC se	RC-SWNIR-MOS Spectrometer head			6	40 to	1050									8	IR-enhanced CMOS linear image sensor

Ultra-compact type	for ir	nstallation into m	obile r	neasui	ing eq	uipme									
Type no.		Туре	200	400	600	800		ral res			2200	2400	2600	Spectral resolution max. (nm)	Image sensor
C10988MA	eries	MS-VIS-MOS Spectrometer head		340	to 750	0								14	CMOS linear
NEW C11708MA	MS se	MS-SWNIR-MOS Spectrometer head			64	40 to	1050							20	image sensor

Information described in this material is current as of March, 2012.

Product specifications are subject to change without prior notice due to improvements or other reasons. Before assembly into final products, please contact us for the delivery specification sheet to check the latest information.

Type numbers of products listed in the delivery specification sheets or supplied as samples may have a suffix "(X)" which means preliminary specifications or a suffix "(Z)" which means developmental specifications.

The product warranty is valid for one year after delivery and is limited to product repair or replacement for defects discovered and reported to us within that one year period. However, even if within the warranty period we accept absolutely no liability for any loss caused by natural disasters or improper product use.

Copying or reprinting the contents described in this material in whole or in part is prohibited without our prior permission.

### MAMATSU

www.hamamatsu.com

HAMAMATSU PHOTONICS K.K., Solid State Division

1126-1 Ichino-cho, Higashi-ku, Hamamatsu City, 435-8558 Japan, Telephone: (81) 53-434-3311, Fax: (81) 53-434-5184 1120-1 ICHINO-CRO, HigdsRin-Ru, Harmamatsu City, 435-858 Japan, Telephone: (81) 53-434-3311, Fax: (81) 53-434-31184
U.S.A.: Hamamatsu Corporation: 360 Foothill Road, P.O.Box 6910, Bridgewater, N.J. 08807-0910, U.S.A., Telephone: (1) 908-231-0960, Fax: (1) 908-231-1218
Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49) 8152-375-0, Fax: (49) 8152-265-8
France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: 33-(1) 69 53 71 00, Fax: 33-(1) 69 53 71 10
United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire AL7 1BW, United Kingdom, Telephone: (44) 1707-294888, Fax: (44) 1707-325777
North Europe: Hamamatsu Photonics Norden AB: Thorshamnsgatan 35 16440 Kista, Sweden, Telephone: (46) 8-509-031-00, Fax: (46) 8-509-031-01
Italy: Hamamatsu Photonics Italia S.R.L.: Strada della Moia, 1 int. 6, 20020 Arese, (Milano), Italy, Telephone: (39) 02-935-81-733, Fax: (39) 02-935-81-741
China: Hamamatsu Photonics (China) Co., Ltd.: 1201 Tower B, Jiaming Center, No.27 Dongsanhuan Beilu, Chaoyang District, Beijing 100020, China, Telephone: (86) 10-6586-6006, Fax: (86) 10-6586-8866