

## DynaSpect Series

DynaSpect is the name for a series of spectrophotometer devices.

Multi-channel Analyzer

# PMA-50



### 【Scientific applications】

UV to visible spectroscopy

Fluorescence spectroscopy

Raman scattering

Chemiluminescence analysis

Liquid chromatography

Gas chromatography

ICP emission analysis

Discharge spectrum analysis

Combustion analysis

Micro spectroscopy

### 【Industrial applications】

Water quality testing

Evaluation of light emitting devices and light sources

Impurities testing

Film thickness measurements

UV radiation measurements

Plasma monitors

Fruit testers

Combustion monitoring

Color filter evaluation

**HAMAMATSU**  
PHOTON IS OUR BUSINESS

# PMA hotonic Multichannel Analyzer

## Simultaneous multi-wavelength measurement !

Measures emission, absorption and reflection spectra in a split second.



The PMA-50 is a multi-channel spectrophotometer consisting of a polychromator, multi-channel photodetector, controller, and data analyzer. The PMA-50 simultaneously measures a wide range of wavelengths without running a wavelength scan, and so can measure emission, absorption and reflection spectra all at one time. Photodetectors are selectable according to measurement wavelengths and light intensity, so you can configure an optimal system that best matches your application.

### FEATURES

#### Simultaneous multi-wavelength measurement

Using a self-scanning image sensor enables simultaneous multi-wavelength measurement without wavelength scanning. Even spectra with weak light emission can be rapidly measured.

#### High-precision spectrum measurement

Ideal for making high-precision spectrum measurements since the axis of the measured spectrum wavelength and photodetector spectral response characteristics can be calibrated.

#### Choice of optimal photodetectors to match application


A wide variety of image sensors are available as photodetectors. Selecting a photodetector that best matches the measurement wavelength range and light intensity makes the PMA-50 the ideal tool for use in a broad range of applications.

#### Up to 3 gratings can be installed at the same time

Different types of gratings are available and up to 3 gratings can be installed at the same time. Even just one polychromator allows you to make simultaneous spectrum observations over a wide range of wavelengths as well as optical measurement with high spectral resolution.


#### Ultra-high sensitivity model BT-CCD system

The BT-CCD system is a multi-channel photodetector using a back-thinned CCD sensor with enhanced sensitivity that makes it ideal for low-light-level measurement. Air-cooled thermoelectric cooling greatly reduces dark current noise so that even weak light can be measured at a high S/N ratio.



#### Near infrared model InGaAs system

These models have an internal InGaAs linear image sensor to allow measuring reflection and absorption spectra in the near infrared region with a wide dynamic range. Two models are available: C10470-256W (256 channels, spectral range 1200 nm to 2600 nm) and C10470-512S (512 channels, spectral range 900 nm to 1700 nm).



**MEASUREMENT MODES**

**Standard measurements**

This is the most basic measurement mode.

Applications: emission spectra for light sources, fluorescence, plasma and the like.

**Reflective measurements**


This is the measurement mode for finding spectral reflectance.

Applications: reflectance measurements for optical filters, coatings and the like.

**Transmittance and absorption measurements**

This is the measurement mode for finding spectral transmittance and absorption.

Applications: measurements of transmittance and absorption in optical filters, films, solutions and the like.



**DISPLAY MODES**

**Spectrum display**

**Display of changes over time**

**3-D display**

**Reflectivity display**

**Transmittance display**

**Absorbance display (OD)**

**Grating list A6751**


No. of grooves (gr/mm)	Blaze wavelength	Spectral range (nm)	Product No.	BT-CCD C10469-1007 (200 nm~1100 nm)				InGaAs C10470-512S (900 nm~1700 nm)				InGaAs C10470-256W (1200 nm~2600 nm)			
				C5094		C5095		C5094		C5095		C5094		C5095	
				Spectral width (nm)	Resolution (nm)	Spectral width (nm)	Resolution (nm)	Spectral width (nm)	Resolution (nm)	Spectral width (nm)	Resolution (nm)	Spectral width (nm)	Resolution (nm)	Spectral width (nm)	Resolution (nm)
50	600 nm	450~1200	40050600	—	—	944	3.2	768	3.9	492	2.5	768	6.0	492	3.84
100	450 nm	250~900	40100450	737	2.5	472	1.6	384	2.0	246	1.2	384	3.0	246	1.92
	780 nm	500~2000	40100780												
150	300 nm	250~800	40150300	492	1.7	315	1.1	256	1.3	164	0.83	256	2.0	164	1.28
	500 nm	300~1800	40150500												
	1250 nm	750~2500	40151250												
300	300 nm	250~1000	40300300	246	0.84	157	0.54	128	0.65	82	0.42	128	1.0	82	0.64
	500 nm	300~1000	40300500												
	750 nm	450~2000	40300750												
	1 μm	650~2400	40301000												
	1.7 μm	1000~2600	40301700												
600	300 nm	250~600	40600300	123	0.42	79	0.27	64	0.33	41	0.21	64	0.5	41	0.32
	500 nm	350~1600	40600500												
	500 nm	360~1250	45060501												
	1000 nm	600~2200	40601000												
830	1500 nm	800~2200	40601500	89	0.30	57	0.19	46	0.23	30	0.15	46	0.36	30	0.23
	1200 nm	650~1600	40831200												
	250 nm	200~550	45120251												
1200	500 nm	350~800	45120501	61	0.21	39	0.13	—	—	—	—	—	—	—	—
	190 nm~700 nm	190~700	45120190												
	1200 nm	600~1100	41201200												
	400 nm~1100 nm	400~1100	45120400												

\*this Spectral range is reference values.


SYSTEM EXAMPLES

Emission spectrum measurements


For fluorescent samples such as fluorescent lamps and EL devices




Fluorescence indicator (Fluorescein) emission spectrum



Chemiluminescence emission spectrum



Fluorescent lamp phosphor emission spectrum



<Configuration>

- Standard PMA-50 configuration

Options


- C10471-01,-02 PMA data analyzer
- Excitation light source: laser, xenon lamp, etc.
- A6751 solution sample holder for transmissivity and fluorescence measurements

<Applications>


- Fluorescence spectroscopy
- Monitoring chemical light emissions

Emission spectrum measurements


Emission spectrum measurements for plasma, electric discharge, ablation and the like




Electric discharge emission spectrum




Emission spectrum during oxide film etching



Emission spectrum during nitride film etching



Temporal changes in plasma emission spectra



<Configuration>

- Standard PMA-50 configuration


Options

- C10471 PMA data analyzer
- A9607 C-mount position checking fiber adapter
- A4869 UV lens

<Applications>

- Plasma component analysis
- Analysis of various emission phenomena

**Absorption spectrum measurements** Spectral transmittance and absorption measurements in optical filters, films, solutions and the like



**<Configuration>**

- Standard PMA-50 configuration


**Options**

- C10471 PMA data analyzer
- L6759 high stability 150 W xenon light source
- A6751 solution sample holder for transmissivity and fluorescence measurements


**<Applications>**

- Absorption spectrum evaluations for solutions and films
- Component analysis for samples
- Monitoring chemical changes


**Component analysis of plastics using transmission spectra (polycarbonate and PET resins)**




**Didymium film absorption spectrum**




**MMA and PMMA transmission spectra**



**Changes in transmissivity in the polymerization from MMA to PMMA (wavelength: 1615 nm)**



**Reflective spectrum measurements** Measurement of spectral reflectance in optical filters, anti-reflective films (AR coatings) and the like




**<Configuration>**

- Standard PMA-50 configuration

**Options**

- C10471 PMA data analyzer
- L6759 high stability 150 W xenon light source

**AR coating reflection spectrum**



**<Applications>**


- Inspection of coatings
- Monitoring thin film growth

## INPUT OPTICS


- Optical probes
  - ① A9080 (BT-CCD)
  - ② A9083 (InGaAs)

Optical probes A9080 and A9083 are optical fibers for guiding light to the polychromator. These optical probes make optical measurement easy because there is no need to set up an optical system. The A9080 and A9083 optical probes are each used with a fiber light condenser A5761-01.

### ① Spectral transmittance of A9080



### ② Spectral transmittance of A9083



- Fiber light condenser A5761-01



The A5761-01 is a fiber light condenser designed to easily attach to the C5094 and C5095 polychrometers. The A5761-01 efficiently inputs measurement light from the optical fiber probe into the polychrometer.

- ▲ A5761-01 connected to an optical fiber probe

## POLYCHROMETERS / POLYCHROMETER ADAPTERS

- Polychromators C5094, C5095



### ▲C5094

The C5094 and C5095 are Czerny-Turner type polychrometers. A variety of gratings are available for each of them and up to 3 gratings can be simultaneously installed. Even one polychromator allows you to make simultaneous spectrum observations over a wide range of wavelengths as well as optical measurement with high spectral resolution.

Since both polychrometers include a GP-IB interface, you can use your PC to switch the gratings, set observation wavelengths, and adjust the entrance slit width.

\* Polychrometers already available at your company can also be used. Please consult our sales office for more information.  
\* Choose a polychromator adapter that matches the detector to be used.

## PHOTODETECTORS

- BT-CCD Multi-Channel detectors C10469-1007




The C10469-1007 is a multi-channel photodetector with a back-thinned CCD sensor that makes it ideal for low-light-level measurement. The C10469-1007 consists of a detector head and a controller. The CCD sensor in the detector head has a high quantum efficiency of 90% or more at a peak wavelength

(approximately 650 nm). Even in the UV region (200 nm to 400 nm), the quantum efficiency is as high as 40% or more. What's more, air-cooled thermoelectric cooling greatly reduces dark current noise for making low-light-level measurements at a high S/N ratio.

Model	C10469-1007
Wavelength	200 nm to 1100 nm
Number of photosensitive device channels	1024 ch
Device cooling temperature	-15 °C
Read-out noise	10 electrons
Dark current	75 electrons/scan (-15 °C:20 ms)

### ● Back-thinned CCD linear image sensor



- InGaAs Multi-Channel detectors C10470-512S
- InGaAs Multi-Channel detectors C10470-256W




The C10470-512S and C10470-256W are high sensitivity multi-channel detectors specifically developed for spectrophotometry in the infrared region. These multi-channel detectors contain an InGaAs image sensor device integrating an InGaAs photodiode array and CMOS signal processing circuit into one package. The C10470-512S uses a single-stage Peltier cooler for thermoelectric cooling and covers a spectral range from 900 nm to 1700 nm.

The C10470-256W, on the other hand, uses a two-stage Peltier cooler and its spectral response ranges from 1200 nm up to 2600 nm.

Model	C10470-512S	C10470-256W
Wavelength	900 nm to 1700 nm	1200 nm to 2600 nm
Number of photosensitive device channels	512 ch	256 ch
Device cooling temperature	-10 °C	-20 °C
Read-out noise	12 500 electrons	12 500 electrons
Dark current	6750 electrons/scan (-10 °C:5 ms)	6.25×10 <sup>6</sup> electrons/scan (-20 °C:5 ms)

### ● InGaAs linear image sensor



**● PMA controller**

The PMA controller drives the C10469-1007, C10470-512S and C10470-256W multi-channel detectors, converts output signals from the sensor into digital signals, and transfers the data to a PC via the USB interface. Measurement conditions are controllable by U6039 PMA basic software that comes with the PMA-50.

**● Control section**

- Interface----- USB 2.0 interface

**● AD converter**

- AD conversion----- 16 bits

**● Functions**

- Data transfer----- Transfers data stored in memory to a PC
- Exposure time setting----- Sets from 5 ms to 64 seconds (depending on detector head)
- Trigger mode setting----- Sets internal sync and external sync
- Trigger polarity setting----- Sets trigger polarity during external sync
- Shutter open/close----- Opens and closes shutter
- Amp gain setting----- Sets data amp gain to 1, 2, or 5 times
- Pixel clock setting----- Sets sensor readout speed

**DATA ANALYZER**



This is an IBM PC compatible computer used for data analysis.

**PMA data analyzer**

- Notebook PC C10471-01
- Desktop PC C10471-02

**PMA SOFTWARE U6039-01**

- Measurement functions ..... Monitoring measurement  
Data measurement
- Temporal resolution measurement functions ... Temporal fluctuation of spectra  
Temporal fluctuation in reflectivity and transmissivity
- Data acquisition condition settings ..... Exposure time settings  
Memory integration count assignment
- Calibration/correction ..... Wavelength axis calibration  
Dark current correction
- Display functions ..... Spectrum display  
Display temporal waveform fluctuations
- Wavelength axis display ..... Wavelength, wave number,  
Raman shift, energy (eV)
- Brightness axis display ..... Linear, logarithm
- Cursor analysis functions ..... Wavelength (wave number, etc.)  
vs. intensity  
Peak detection  
FWHM measurement  
Integrated intensity
- Other analytical functions ..... Smoothing  
Differential waveform

**OPTIONS**



**Solution sample holder for transmissivity and fluorescence measurements A6751**

This is a dedicated holder with an integrated condensing lens for the use with vials.



**Reflection measurement optics A9665**

These are optics making it possible to illuminate the sample at 45° from the light source and measure the reflected light.



**C-mount fiber adapter A6399**

This is an adapter for securing the fiber input optics to the C-mount of a microscope or the like.



**C-mount position checking fiber adapter A9607**

In addition to the function of the C-mount fiber adapter, the measurement position can be checked. Measurements in the UV.




**UV lens A4869**

Condensing lens for UV. f=50 mm, F3.5 (A6399, A8482 required)




**Light source L5351 for wavelength calibration**




**High stability 150 W xenon light source L6759**

This is a high stability xenon light source with output wavelengths from 250 nm to 1000 nm for excitation and absorption measurements.



**Halogen light source L6758**

This is a halogen light source with output wavelengths from 400 nm to 1000 nm for excitation and absorption measurements.



SPECIFICATIONS


	PMA-50/BT-CCD system	PMA-50/InGaAs system	
<b>Model</b>	Ultra-high sensitivity	Near infrared model	
Photo-detector	BT- CCD linear image sensor C10469-1007	InGaAs linear image sensor C10470-512S	InGaAs linear image sensor C10470-256W
Wavelength	200 nm to 1100 nm	900 nm to 1700 nm	1200 nm to 2600 nm
Number of photosensitive device channels	1024 ch	512 ch	256 ch
Device cooling temperature	-15 °C	-10 °C	-20 °C
Read-out noise	10 electrons	12 500 electrons	12 500 electrons
Dark current	75 electrons/scan (-15 °C : 20 ms)	6750 electrons/scan (-10 °C : 5 ms)	6.25 × 106 electrons/scan (-20 °C : 5 ms)
AD resolution	16 bit		

● Polychrometers


Model	C5094	C5095
Optical layout	Czerny-Turner type (comes with aberration-corrected toroidal mirror)	
Focal length	250 mm	500 mm
F value	4	8
Entrance slit width	Variable from 10 μm to 2000 μm	
Grating	Up to 3 gratings can be installed at the same time.	
Reciprocal linear dispersion	2.5 nm/mm (at 1200 gr/mm)	1.6 nm/mm (at 1200 gr/mm)

**DIMENSIONAL OUTLINES (Unit :mm)**


● Fiber input optics (approx.100 g)




● Polychrometers C5094 (approx.13.5 kg)




● Polychrometers C5095 (approx.16 kg)




● Controller (approx.3 kg)




● Detector C10469-1007 (approx. 520 g)



● Detector C10470-512S (approx. 940 g)



● Detector C10470-256W (approx. 870 g)



- ★ DynaSpect is a registered trademark of Hamamatsu Photonics K.K.
- ★ Product and software package names noted in this documentation are trademarks or registered trademarks of their respective manufacturers.
- ★ The spectral response graphs in this document are representative examples, so they are not guaranteed.
- Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications and external appearance are subject to change without notice.

© 2011 HAMAMATSU PHOTONICS K.K.

**HAMAMATSU**

Web site [www.hamamatsu.com](http://www.hamamatsu.com)

HAMAMATSU PHOTONICS K.K., Systems Division

812 Joko-cho, Higashi-ku, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-435-1574, E-mail: export@sys.hpk.co.jp

U.S.A. and Canada: Hamamatsu Corporation: 360 Foothill Road, Bridgewater, N.J. 08807-0910, U.S.A., Telephone: (1) 908-231-0960, Fax: (1) 908-231-0852, E-mail: usa@hamamatsu.com

Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-2658, E-mail: info@hamamatsu.de

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10, E-mail: infos@hamamatsu.fr

United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road, Welwyn Garden City, Hertfordshire, AL7 1BW, U.K., Telephone: (44) 1707-294888, Fax: (44) 1707-325777, E-mail: info@hamamatsu.co.uk

North Europe: Hamamatsu Photonics Norden AB: Smidesvägen 12, SE-171-41 Soina, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01, E-mail: info@hamamatsu.se

Italy: Hamamatsu Photonics Italia S.R.L.: Strada della Moia, 1 int.6-20020 Arese (Milano), Italy, Telephone: (39)02-935 81 733, Fax: (39)02-935 81 741, E-mail: info@hamamatsu.it

Cat. No. SDSS 0003E02

AUG/2011 HPK

Printed in Japan